Interpreter Version History

2.26	
	Fixed copy bugs, DEMO segment export disabled, change license on Reviewer and Interpreter DEMO
2.25	
	Fixed segment export bugs
2.24_	<u></u>
	Fixed debug exception bug upon cancellation of segment export Fixed demo annotation bug from PC-Display
2.23	
	Fixed Celsius display bug
-	Added Segment export feature
2.22_	<u> </u>
_	Fixed Interpreter bug which created artifact in CO2 waveform when
-	adding new beat markers Fixed Go to Beat command to correctly find the first beat (Beat 1) in long record where Beat 1 is not at the start of the file.
2.21	
	Fixed Reviewer bug which showed QuickChek controls when changing from CO2 back to ECG waveform review
_	Updated Interpreter manual
-	Darkened CO2 traces
_	Fixed deployment bug with Species&Breeds file
2.20	
	Fixed debug assertion error for QuikChek
-	Changed CO2 waveform width and color
0.10	
2.19_	Added CO2 waveform display feature
_	
_	
_	Added CO2 scale pull down and default setting
-	Added watermarks for demo and Tilley episodes
2.16	
	Tixed index into episode clicking on an event bug
-	
2.15_	<u></u>
-	Fixed E-mail Icon bug
_	Added pull down menu and customizable patient species and breed
2.14_	added persional consimumation leadersto
_	Added Banfield configuration lockouts Added HTML manual and link to support website
	Fixed registration bug
2.13_	<u></u>
	Migrated registration key to HKLM

а

2		12	
	•	12	

- Changed method of opening Tilley example episodes

___2.11___

- New 24 bit color icons on the toolbar
- Report formatting change to use landscape printing for Vital signs to reduce the number of pages in the printout
- Added "Vitals Only" button to the toolbar to reduce the number of events displayed
- Added an Export to Excel
- Added PC-Display style logo to reports.
- Added International clinic Vet Info entry
- Updated manual
- Ported to Visual C 2005 environment.
- Switched to MSI style installer.